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[Osteoblast response to continuous phase macroporous scaffolds under static and dynamic culture conditions.](#) Meretoja VV, Malin M, Seppälä JV, Närhi TO. *J Biomed Mater Res.* 2008;89A(2):317-25.

[Molecular basis for action of bioactive glasses as bone graft substitute.](#) Välimäki VV, Aro HT. *Scand J Surg.* 2006;95(2):95-102.

[Intact surface of bioactive glass S53P4 is resistant to osteoclastic activity.](#) Wilson T, Parikka V, Holmbom J, Ylänen H, Penttinen R. *J Biomed Mater Res.* 2005;77A(1):67-74.

[Granule size and composition of bioactive glasses affect osteoconduction in rabbit.](#) Lindfors NC, Aho AJ. *J Mater Sci: Mater Med.* 2003;14(4):265-372.

[Osteoblast differentiation of bone marrow stromal cells cultured on silica gel and sol-gel-derived titania.](#) Dieudonné SC, van den Dolder J, de Ruijter JE, Paldan H, Peltola T, van 't Hof MA, Happonen RP, Jansen JA. *Biomaterials.* 2002;23(14):3041-51.

[Histomorphometric and molecular biologic comparison of bioactive glass granules and autogenous bone grafts in augmentation of bone defect healing.](#) Virolainen P, Heikkilä J, Yli-Urpo A, Vuorio E, Aro HT. *J Biomed Mater Res.* 1997;35(1):9-17.

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[In vitro antibacterial activity of bioactive glass S53P4 on multiresistant pathogens causing osteomyelitis and prosthetic joint infection.](#) Cunha TM, Murça MA, Nigro S, Klautau GB, Salles MJC. *BMC Infectious Diseases.* 2018;18:157.

[Recent evidence on bioactive glass antimicrobial and antibiofilm activity: a mini-review.](#) Drago L, Toscano M, Bottagisio M. *Materials.* 2018;11(2):326.

[Influence of bioactive glass S53P4 granules and putty on osteomyelitis associated bacteria in vitro.](#) Stoor P, Frantzen F. *Biomed Glasses.* 2017;3(1).

[Antibiofilm agents against MDR bacterial strains: Is bioactive glass BAG-S53P4 also effective?](#) Bortolin M, De Vecchi E, Romanò CL, Toscano M, Mattina R and Drago L. *J Antimicrob Chemother.* 2016 Jan;71(1):123-7.

[Bioactive glass combined with bisphosphonate provides protection against biofilms formed by the periodontal pathogen *Aggregatibacter actinomycetemcomitans*.](#) Hiltunen AK, Skogman ME, Rosenqvist K, Juvone H, Ihalainen P, Peltonen J, Juppola A, Fallalero A. *Int J Pharm* 2016;501:211-20.

[The antibacterial activity of bioactive glass.](#) Trishala A, Jacob C. *International Journal of Advanced Research.* 2016;4(6):1070-7.

[Antimicrobial activity and resistance selection of different bioglass S53P4 formulations against multidrug resistant strains.](#) Drago L, De Vecchi E, Bortolin M, Toscano M, Mattina R and Romanò CL. *Future Microbiol.* 2015;10(8):1293-9.

[Efficacy of antibacterial bioactive glass S53P4 against *S. aureus* biofilms grown on titanium discs in vitro.](#) Coraça-Huber DC, Fille M, Hausdorfer J, Putzer D, Noglén M. *J Orthop Res.* 2014;32(1):175-7.

[In vitro antibiofilm activity of bioactive glass S53P4.](#) Drago L, Vassena C, Fenu S, De Vecchi E, Signori V, De Francesco R, Romanò CL. *Future Microbiol.* 2014;9(5):593-601.

[The in vitro antibacterial effect of S53P4 bioactive glass and gentamicin impregnated polymethylmethacrylate beads.](#) Gergely I, Zazgyva A, Man A, Zuh SG, Pop TS. *Acta Microbiol Immunol Hung.* 2014;61(2):145-60.

[Antibacterial effects and dissolution behavior of six bioactive glasses.](#) Zhang D, Leppäranta O, Munukka E, Ylänen H, Viljanen MK, Eerola E, Hupa M, Hupa L. *J Biomed Mater Res.* 2010;93A(2):475-83.

[Bactericidal effects of bioactive glasses on clinically important aerobic bacteria.](#) Munukka E, Leppäranta O, Korkeamäki M, Vaahtio M, Peltola T, Zhang D, Hupa L, Ylänen H, Salonen JI, Viljanen MK, Eerola E. *J Mater Sci: Mater Med.* 2008;19(1):27-32.

[Antibacterial effect of bioactive glasses on clinically important anaerobic bacteria in vitro.](#) Leppäranta O, Vaahtio M, Peltola T, Zhang D, Hupa L, Ylänen H, Salonen JI, Viljanen MK, Eerola E. *J Mater Sci: Mater Med.* 2008;19(2):547-51.

[Factors controlling antibacterial properties of bioactive glasses.](#) Zhang D, Munukka E, Hupa L, Ylänen H, Viljanen MK, Hupa M. *Key Eng Mat.* 2007;330-332:173-6.

[Effect of ion release on antibacterial activity of melt-derived and sol-gel-derived reactive ceramics.](#) Vaahtio M, Munukka E, Leppäranta O, Zhang D, Eerola E, Ylänen H, Peltola T. *Key Eng Mat.* 2006;309-311:349-52.

[Comparison of antibacterial effect on three bioactive glasses.](#) Zhang D, Munukka E, Leppäranta O, Hupa L, Ylänen H, Salonen J, et al. *Key Eng Mat.* 2006;309-311:345-8.

- [A comparative study on the disinfection potentials of bioactive glass S53P4 and calcium hydroxide in contra-lateral human premolars ex vivo.](#) Zehnder M, Luder H, Schätzle M, Kerosuo E, Waltimo T. International Endodontic J. 2006;39(12):952-8.
- [Interactions between the bioactive glass S53P4 and the atrophic rhinitis-associated microorganism Klebsiella ozaenae.](#) Stoor P, Söderling E, Grenman R. J Biomed Mater Res. 1999;48:869-74.
- [Antibacterial effects of a bioactive glass paste on oral micro-organisms.](#) Stoor P, Söderling E, Salonen JI. Acta Odontol Scand. 1998;56:161-5.
- [Interactions between BAG and periodontal pathogens.](#) Stoor P, Kirstilä V, Söderling E, Kangasniemi I, Herbst K, Yli-Urpo A. Microbial Ecology in Health and Disease. 1996;9:109-14.
- [Interactions between the frontal sinusitis-associated pathogen Haemophilus Influenzae and the bioactive glass S53P4.](#) Stoor P, Söderling E, Andersson OH, Yli-Urpo A. Bioceramics. 1995;8:253-8.

SEPTISCHE KNOCHENCHIRURGIE

- [Clinical applications of Bioactive glass S53P4 in bone infections: a systematic review.](#) Bigoni M, Turati M, Zanchi N, Lombardo AS, Graci J, Omeljaniuk RJ, Zatti G, Gaddi D. Eur Rev Med Pharmacol Sci. 2019;23(2 Suppl):240-51. [Open access]
- [The use of bioactive glass S53P4 as bone graft substitute in the treatment of chronic osteomyelitis and infected non-unions – a retrospective study of 50 patients.](#) Al Malat TA, Glombitza M, Dahmen J, Hax P-M, Steinhausen E. Z Orthop Unfall 2018;156:152-9.
- [Bioactive glass S53P4: a new opportunity for the treatment in the diabetic foot osteomyelitis.](#) De Giglio R, Stefani I, Mondello T, De Filippis G, Mazzone A. Eur J Intern Med. 2018 Aug;54:e15-e16.
- [Contemporary treatment of chronic osteomyelitis: implementation in low- and middle-income countries.](#) Geurts J, Vranken T, Gabriels F, Arts JJ, Moh P. SA Orthop J. 2018;17(2):40-3.
- [Evaluation of the clinical effectiveness of bioactive glass \(S53P4\) in the treatment of non-unions of the tibia and femur: study protocol of a randomized controlled non-inferiority trial.](#) Tanner MC, Heller R, Westhauser F, Miska M, Ferbert T, Fischer C, Gantz S, Schmidmaier G, Haubruck P. Trials. 2018;19(1):299.
- [Treatment of cavitary bone defects in chronic osteomyelitis: Biogactive glass S53P4 vs. calcium sulphate antibiotic beads.](#) Ferrando A, Part J, Baeza J. Journal of Bone and Joint Infection 2017;2(4):194-201.
- [Il trattamento delle pseudoartrosi settiche di tibia con la tecnica di Ilizarov isolata o associata a biovetro attivo. Treatment of infected tibial nonunions with Ilizarov technique and bioactive glass.](#) Testa G, Aloj DC, Pavone V, Massè A, Avondo S, Sessa G. Giornale Italiano di Ortopedia e Traumatologia 2017;43:251-7.
- [Antibacterial bioactive glass, S53P4, for chronic bone infections - a multinational study.](#) Lindfors NC, Geurts J, Drago L, Arts JJ, Juutilainen V, Hyvönen P, Suda A, Domenico A, Artiaco S, Alizadeh C, Brychcy A, Bialecki J, Romano C. Adv Exp Med Biol. January 2017;971:81-92.
- [Reconstruction of vertebral bone defects using an expandable replacement device and bioactive glass S53P4 in the treatment of vertebral osteomyelitis: three patients and three pathogens.](#) Kankare J, Lindfors NC. Scand J Surg. 2016 Feb 29.
- [Two-stage treatment for infected nonunion of the humerus: a case report.](#) Puma Pagliarello C, Testa G, Lombardini A, Condorelli G, Pavone V. Infec Dis Trop Med 2016;3(2):E390.
- [Clinical applications of S53P4 bioactive glass in bone healing and osteomyelitic treatment: a literature review.](#) van Gestel NA, Geurts J, Hulsen DJ, van Rietbergen B, Hofmann S, Arts JJ. Biomed Res Int. 2015;2015: 684826.
- [Clinical application of antimicrobial bone graft substitute in osteomyelitis treatment: a systematic review of different bone graft substitutes available in clinical treatment of osteomyelitis.](#) van Vugt TA, Geurts J, Arts JJ. Biomed Res Int. 2016; Article ID 6984656.
- [Treatment of osteomyelitis by means of bioactive glass - initial experience in the Netherlands.](#) Geurts J, Vranken T, Arts JJ. Nederlands Tijdschrift voor Orthopaedie. 2016;23(2):37-41.
- [Bioactive glass for long bone infection: a systematic review.](#) Aurégan J-C, Bégué T. Injury, Int. J. Care Injured. 2015;46(8):3-7.
- [New trends in surgical treatment of chronic osteomyelitis.](#) Juutilainen V, Lindfors NC. Suomen Ortopedia ja Traumatologia. 2015;38:59-61.
- [A comparative study of the use of bioactive glass S53P4 and antibiotic-loaded calcium-based bone substitutes in the treatment of chronic osteomyelitis - a retrospective comparative study.](#) Romanò CL, Logoluso N, Meani E, Romanò D, De Vecchi E, Vassena C, Drago L. Bone Joint J 2014;96-B:845-850.
- [Bioactive glass BAG-S53P4 for the adjunctive treatment of chronic osteomyelitis of the long bones: an in vitro and prospective clinical study.](#) Drago L, Romanò D, De Vecchi E, Vassena C, Logoluso N, Mattina R, Romanò CL. BMC Infectious Diseases 2013;13:584.
- [Review: Emerging developments in the use of bioactive glasses for treating infected prosthetic joints.](#) Rahaman M, Bal BS, Huang W. Mater Sci Eng C Mater Biol Appl. 2014;C41:224-31.
- [Through the looking glass; bioactive glass S53P4 \(BonAlive®\) in the treatment of chronic osteomyelitis.](#) McAndrew J, Efrimescu C, Sheehan E, Niall D. Ir J Med Sci. 2013;182(3):509-511.
- [Posttraumatic osteomyelitis.](#) Juutilainen V. Suomen Ortopedia ja Traumatologia 2011;34:38-41.
- [Bioactive glass S53P4 as bone graft substitute in treatment of osteomyelitis.](#) Lindfors NC, Hyvönen P, Nyyssönen M, Kirjavainen M, Kankare J, Gullichsen E, Salo J. Bone. 2010;47:212-218.

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Bonalive Biomaterials Ltd